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
**Lanzetten** STERILE R

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**Manufacturer:**  
**Shandong Lianfa Medical Plastic Products Co.,Ltd.**  
 No.1 Shuangshan Sanjian Road, 250200 Zhangqiu City, Jinan, Shandong, P.R. China  
 Tel: +86 531 6132-8777    linkfar@sdlianfa.com

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**EPS BIO Technology Corp.**  
**No.8, R&D RD. III, Hsinchu Science Park, Hsinchu, Taiwan 30076**

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





- 1 Display
- 2 Battery compartment (on the back)
- 3 > - button for memory access, to enter values and to look at test results
- 4 ⏻ - button, to confirm settings, delete results or to switch on the device (press and hold for approx. 3 seconds).
- 5 Insertion slot for test strips
- 6 Test strip
- 7 Connection for USB cable

- 8 Blood sample area (absorbent slot)
- 9 Reaction cell
- 10 Grip/holder
- 11 Insert in the direction of the arrow
- 12 Contact electrode
- 13 AST cap
- 14 Sleeve
- 15 Hub
- 16 Trigger button
- 17 Cap end (adjustable)

- 19 Time (hours:minutes)  
AM = before noon; PM = after noon
  - 20 Memory slot
  - 21 Blood glucose reading
  - 22 Unit of measure
  - 23 Battery symbol (low battery)
  - 24 Ambient temperature error
  - 25 Insert test strip
  - 26 Apply blood or control solution
  - 27 Active USB connection
  - 28 System error
  - 29 Warning regarding a possible diabetical ketoacidosis.  
Consult doctor!
  - 30 Control solution measurement mode
  - 31 Date (month:day)
  - 32 Alarm (before meals with automatic alarm after 2 hours to remind you to perform a test after meals)
- Before meals (AC)
- After meals (PC)

The following signs and symbols on items, packaging and in the instruction manual bear important information:

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

#### WARNING

These warning notes must be observed to prevent any injury to the user.

#### CAUTION

These notes must be observed to prevent any damage to the device.

#### NOTE

These notes give you useful additional information on the installation or operation.

LOT

Lot number

REF

Catalogue number

SN

Serial number

Manufacturer

Recycling symbols/codes:

These are used to provide information about the material and its proper use and recycling.

Ambient pressure limitation

Humidity range

IVD In vitro diagnostic medical device (for external use only)

The blood glucose meter is in compliance with Directive 98/79/EC for in vitro diagnostic devices.

Do not re-use

STERILE R

Radiation sterilized

Use by

Biological risks

MD

Medical Device

Temperature range

Date of manufacture

YYYY-MM-DD

Contains sufficient for <n> tests

CONTROL

Control solution

Discard 6 months after opening

3V (CR2032) x 2-battery only

Keep away from sunlight

EC REP

Authorized representative in the European Community

mg/dL

Blood glucose test result in mg/dL

mmol/L

Blood glucose test result in mmol/L

Large LCD screen

5-seconds result

Hypo alarm

0.6 µL blood volume

Alternative site testing

Human-factor flow

Importer

### IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!

**Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.**

**Should you give this device to another person, it is vital that you also pass on these instructions for use.**

#### Intended Use

The MediTouch 2 Self Monitoring Blood Glucose Test System is intended for the quantitative measurement of glucose in fresh capillary whole blood from fingertips, palm and forearm. Testing is done outside the body (In Vitro diagnostic use). It is indicated for self-testing (over the counter [OTC]) by persons with diabetes, or in clinical settings by healthcare professionals, as an aid to monitor the effectiveness of diabetes control.

#### Contraindications

The system is not suitable for diagnosing diabetes for children younger than 12 years. For use on older children ask your doctor.

It is not suitable for diagnosing diabetes or for testing the blood glucose levels of newborn babies.

#### Always observe the following points

- Although the MediTouch 2 SMBG System is easy to use, you may need to consult with your healthcare professional (this may be your doctor, pharmacist or diabetes nurse educator) for instructions on how to use the system. Only the correct use of the system will ensure accurate results.
- This device measures in mg/dL.
- The unit is only to be used for the specific purpose described in this instruction manual.
- Any misuse will void the warranty.
- MediTouch 2 Blood Glucose meter is designed and approved for testing fresh capillary blood from fingertip, palm and forearm. The meter is for in vitro diagnostic use ONLY (for testing outside the body). It should not be used to diagnose or screen for diabetes.
- MediTouch 2 Blood Glucose Meter can only be used with MediTouch 2 Blood Glucose Test Strips. Other test strips will give inaccurate results.
- Only use the accessories which have been recommended by the manufacturer (test strips, control solution) with this system.
- This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device.
- Children must be supervised to ensure that they do not play with the device.
- Do not operate the device in the vicinity of high-frequency transmitters, e.g. microwave and short wave transmitters.
- Do not use the device if it is not working correctly, if it has been dropped or has fallen into water or has been damaged.
- Protect the unit from moisture. Should moisture enter the unit, remove the batteries and stop using it immediately. Contact your authorised service centre.
- If a fault occurs, do not try to repair the unit yourself. Attempts to do so will void the warranty. Refer all servicing to authorized service personnel.
- Always keep the meter clean and store it in a safe place. Protect the meter from direct sunlight to prolong its service life.
- Do not store the meter and the test strips in a vehicle, bathroom or refrigerator.
- Extremely high humidity can affect the test results. A relative humidity of more than 90% can cause incorrect readings.
- Store the meter, the test strips and the lancing device out of reach of children and pets.
- Remove the battery if you do not intend to use the meter for one month or longer.
- Consult your healthcare professional before making a blood glucose test with this device.
- Warning for potential biohazard: Healthcare professionals using this system on multiple patients should be aware that all products or objects that come in contact with human blood, even after cleaning, should be handled as if capable of transmitting a viral disease.

- Consult with your healthcare professional before testing on your palm or forearm.
- Do not touch the strips with wet hands.
- Do not use expired strips (the expiration date is shown on the bottle.)
- Do not bend, cut or twist the strips.
- Altitude up to 3,048 meters above sea level has no effect on readings.

#### Instructions for your health

- This blood glucose meter is intended for actively testing a person's blood glucose at home. The reading from a blood glucose home-test system does not replace a professional test performed in a laboratory.
- The meter is intended for external use only (in vitro).
- Only use fresh, capillary whole blood from the finger tips for the test.
- You may only adapt the procedure for using products at home and self-monitoring, if you have first received the appropriate training to do so.
- Do not change any treatments as a result of your blood glucose reading without first consulting your doctor.
- The system is not suitable for testing critically ill patients.
- Your meter only requires one small drop of blood to perform a test. You can get this from a fingertip. Use different place for each test. Repeated lancing in the same place can lead to infection and numbness.
- Test results below 60 mg/dL (3.3 mmol/L) are an indication of hypoglycemia, meaning the blood glucose level is too low. If the reading is above 240 mg/dL (13.3 mmol/L), symptoms of a high blood glucose level (hyperglycemia) can occur. Consult your doctor, if your readings are regularly above or below these levels.
- If the test results display "HI" or "LO", perform the test once again. If you obtain similar readings again, consult your doctor immediately and follow his instructions.
- If the proportion of red blood cells (haematocrit level) is very high (over 60 %) or very low (under 20 %), this could distort the test results.
- Dehydration or a lack of water (such as from sweating) can result in incorrect readings. If you think you are suffering from dehydration or a lack of water, consult your doctor immediately.
- If you have followed the instructions in this manual and symptoms persist which are not associated with your blood glucose level or your blood pressure, consult your doctor.
- For additional advice on your health, read the instruction manual for the test strips carefully.

#### WARNING - Risk of infection

- Used test strips and lancets are considered to be hazardous, biological non-biodegradable waste. Carelessness when disposing of these items can lead to the spread of infection. If necessary, consult your local waste disposal company, your doctor or pharmacist.
- Dispose of your used test strips and used lancets carefully. If you dispose of the used parts with the household rubbish, make sure you wrap them first to avoid harming or infecting other people.
- Medical staff or others who use this meter on more than one patient must be aware that all products or objects which come into contact with human blood must be dealt with as if they are capable of spreading viral diseases, even after cleaning.
- Never use a lancet or the lancing device on more than one person.
- Use a new sterile lancet and a new test strip for every test.
- Lancets and test strips are disposable.
- Avoid getting hand cream, oil or dirt in or on the lancet, lancing device and test strip.

- If you are experiencing dehydration, frequent urination, low blood pressure, shock or hyperosmolar hyperglycemic nonketotic coma (HHNC), you may get a test result that is lower than what your blood glucose really is. If you think you are dehydrated, call your doctor right away.
- Please read your test strip instructions carefully for additional health-related information.

#### Warning for potential biohazard

**Healthcare professionals using this system on multiple patients should handle all products or objects in contact with human blood carefully to avoid transmitting viral disease, even after cleaning.**

#### Instructions for using the blood glucose test strips

- Only use these with the MediTouch 2 meter.
- Store the test strips in their original container.
- To avoid contamination, only touch the test strips with clean, dry hands. When removing the strips from the container and inserting them in the meter, only touch them using the grip (holder).
- Close the container again immediately once you have removed the test strip. This keeps the test strips dry and free of dust.
- Use the test strip within three minutes of removing it from the container.
- The test strip is intended to be used only once. Do not use it again.
- Write the opening date on the label of the container when you open it for the first time. Observe the expiry date. The test strips can be used for approx. six months after opening the container or until the expiry date, whichever comes first.
- Run a control solution test every time you open a new box of test strips.
- Do not use test strips which have already exceeded the expiry date, as this can distort the test result. The expiry date is printed on the container.
- Store the test strips in a cool, dry place but not in the refrigerator.
- Store the test strips between 2°C and 30°C (35.6°F - 86°F). Do not freeze the test strips.
- Protect the strips from damp and direct sunlight.
- Do not apply the blood sample or control solution to the test strip before you have inserted it in the meter.
- Only apply the blood sample or the control solution provided on the test strip slot. Applying any other substance will lead to an imprecise or incorrect reading.
- The test strips can be used at altitudes of up to 3,048 m without having any impact on test results.
- Do not bend, cut or adjust the test strips in any way.
- Keep the container with the test strips away from children. There is a risk of choking from the lid. The lid also contains desiccative substances which could be dangerous if they are inhaled or swallowed. This could also lead to skin and eye irritations.

#### Instructions on how to use the control solution

- Only use MediTouch® 2 control solution.
- Only use with the MediTouch® 2 test strips.
- Write the opening date on the label of the container. The control solution can be used for approx. three months after opening the container or until the expiry date, whichever comes first.
- Do not use the control solution after the expiry date.
- The ambient temperature for using the control solution should lie between 10°C - 40°C (50°F - 104°F).
- The temperatures for storing and transporting the control solution should lie between min 2°C and max 30°C (35.6°F - 86°F). Do not store the solution in the refrigerator and do not freeze it.
- Shake the bottle with the control solution well before you open it. Wipe away the first drop and use the second one to ensure a good sample for a precise test result.
- To ensure the control solution does not get contaminated, wipe away the remaining solution on the tip of the container with a clean cloth before closing it again.
- The control solution may stain your clothing. Rinse the soiled clothing with water and detergent.
- Do not put any excess control solution back in the container.
- Close the container carefully after every use.

#### Items supplied and packaging

Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use and contact your supplier or your service centre.

The following parts are included:

- 1 MediTouch® 2 Blood Glucose Meter
- 1 AST Lancing Device Cap
- 1 Lancing device
- 2 CR2032 Lithium Batteries
- 10 MediTouch® 2 Blood Glucose Test Strips
- 1 Self-Test Log Book
- 10 Lancets
- 1 Storage bag
- 1 MediTouch® 2 Level 2 Control Solution 4 ml
- 1 Instruction manual
- 1 Teststrips manual
- 1 Control solution manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your supplier without delay.

**WARNING**  
Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!

- Optional Accessories:
- Level 1 Control Solution
  - Level 3 Control Solution

- Level 2 control solution is included with the system.
- Level 1 and Level 3 control solutions are available on the market. For purchase, please contact your local dealer.

### Special features of the MediTouch® 2 blood glucose meter


Measuring your blood glucose level regularly can be a great help when dealing with your diabetes. This blood glucose meter is designed so that you can use it easily, regularly and anywhere you choose. The lancing device can be set to the sensitivity of your skin. The diabetic diary provided helps you to recognise and record the affect of your eating habits, sport activities or medication on the test results. Always consult your doctor about your test results and treatment. This meter is intended for actively testing a person's blood glucose at home. It is not suitable for diagnosing diabetes or for testing the blood glucose of newborn babies.

Your MediTouch® 2 blood glucose meter consists of five main parts: the blood glucose meter, the lancing device, the lancets, the test strips and the control solution.

These parts are specially designed to be used together and for their quality to ensure precise test results. Only use MediTouch® 2 approved test strips, lancets and control solution for your blood glucose meter. Precise test results can only be ensured when the meter is used properly. Only use fresh, capillary whole blood for the test, preferably from the finger tips.

The meter measures the blood glucose level very precisely. It has an automatic memory for 480 readings with the date and time. It also calculates the average level using the blood glucose readings from the previous 7, 14, 30 and 90 days. This allows you to monitor any changes and discuss them with your doctor.

### Inserting the batteries


1. Open the battery compartment on the back side of the device by sliding the lid in the direction of the arrow.
2. Insert two batteries as shown. The device will emit a short beep after correct insertion.
3. Close the lid of the battery compartment. It will "snap" into its place.

### WARNING BATTERY SAFETY INFORMATION

- Do not disassemble batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Only replace with batteries of the same type, never use different types of batteries together or used batteries with new ones!
- Insert the batteries correctly, observing the polarity!
- Remove the batteries from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge these batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!
- Please observe our disposal instructions in the chapter 'Disposal instructions'!

### Setting the time and date

Setting the current time and date in your meter is important if you intend to use the meter memory.

1. After insertion of the batteries the meter switches on automatically.
2. The last two digits of the year flash in the display. Press > to set the correct value and press to confirm.
3. Repeat step 2 to enter the date and time. The area in which a setting is currently being made flashes.
4. The icon appears in the display. The device is now ready for the test with control solution.

The unit of measure (mg/dL or mmol/L) is installed in the device. If you want to change the unit of measure, contact the customer service.

### Using the control solution

#### Why do I need to perform a test with control solution?

When you perform a test with control solution, you will find out whether your meter and test strips work properly and deliver exact readings. You should perform a test in the following cases:

- You are using the MediTouch 2 Blood Glucose Meter for the first time.
- You open a new bottle of test strips.
- You think the meter or test strips may be working incorrectly.
- You dropped the meter.
- You have repeated a test and the test results are still lower or higher than expected.
- You are practicing the test procedure.

**Professional users are instructed to follow federal, state and local guidelines.**


### Performing a test with control solution

**WARNING** Before you perform a test with control solution, be sure to read the safety instructions in full, in particular the sections 'Instructions for using the blood glucose test strips' and 'Instructions for using the control solution'.


### Operating

You need the device, a test strip and the control solution.

**The measurement results in the control solution measurement mode will not be recorded in the memory for blood glucose readings.**

1. Insert a test strip into the meter in the direction of the arrow. The symbol appears automatically.
2. Press the > - button, until the symbol (control solution measurement mode) appears. The drop symbol flashes. Press to confirm the setting.
3. Place the device on a flat and even surface, e.g. like a table.
4. Remove the cap of the control solution bottle and wipe the tip of the bottle with a tissue.

5. Press the container so that a tiny droplet forms on the tip of the container.
6. Place the drop onto the blood sample area on the end of the test strip.
7. Do not get any control solution on the top of the test strip.
8. When a sufficient amount of control solution has been absorbed by the reaction cell you will hear a beep and in the display "OK" will be shown.
9. The meter starts a countdown of approx. 5 seconds which is shown in the display.
10. In the display, a test result appears. Before you remove the test strip, check whether the result lies within the range indicated on the test strips container.
11. After that, remove the test strip and throw it away.

### Evaluating the control test result

The permissible range for the control solution reading is indicated on the label of the test strip container. Your test result must lie within the range indicated. Make sure that you compare the test result with the correct range. If the control test result lies within the range indicated on the test strip container, then the meter and the test strips are working accurately.

If the control test result does not lie within the range indicated on the test strip container, then the following options will be displayed to rectify the problem:

**Cause:** Has the test strip been lying around open for a long time?  
**Remedy:** Repeat the test with a test strip that has been stored correctly.

**Cause:** Was the test strip container properly closed?  
**Remedy:** The test strips are damp. Replace the test strip.

**Cause:** Did the meter work properly?  
**Remedy:** Repeat the test as described. If problems persist, get in touch with the service centre.

**Cause:** Is the control solution soiled or has it exceeded the expiry date?  
**Remedy:** Use new control solution to check the performance of the meter.

**Cause:** Were the test strips and control solution stored in a cool, dry place?  
**Remedy:** Repeat the control test using strips and solution that have been stored correctly.

**Cause:** Did you follow the steps of the test procedure properly?

**Remedy:** Repeat the test as described. If problems persist, get in touch with the service centre.

### Preparing the blood glucose test

#### Using the lancing device

The lancing device enables you to hygienically and easily draw a drop of blood for the blood glucose test and it is quick and painless. The lancing device can be set to the sensitivity of your skin. You can adjust the tip to 5 different lancing depths. Twist the cap end in the appropriate direction until the arrow is pointing to the number for the lancing depth you want.


- 1) Consult the following for the suitable lancing depth:  
1 - 2 for soft or thin skin, 3 - 4 for normal skin and 5 for thick or callous skin.
- 2) Never use a lancet or a lancing device that belongs to another person.  
Use a new sterile lancet for each test.

**WARNING** Lancets are intended to be used only once. Used test strips and lancets are considered hazardous, biological non-biodegradable waste. Dispose of them taking into account that they are capable of spreading infection. Dispose of the lancets so that there is no risk of injury or infection to other people.

### Inserting a lancet into the lancing device


Before using the lancing device you need to insert a lancet.

**WARNING** Before carrying out a blood glucose test and before using the lancing device, it is essential that you read the safety instructions in full, in particular the sections 'Information for your health' and 'Instructions for using the blood glucose test strips'.

1. Wash your hands with soap and warm water. Rinse and dry thoroughly. If needed, you can also rub the area of skin you have chosen for the blood sample with a special cleansing pad.
2. Open the lancing device by twisting the protective cap in a clockwise direction and then remove it. Insert the lancet all the way (without turning it) into the lancing device. Twist off the protective cap on the lancet carefully.
3. Replace the protective cap on the lancing device and tighten by turning it in an anticlockwise direction. To use the AST cap, check chapter "Using the AST cap".

Markings for lancing depths

4. Set the appropriate lancing depth, as described before.
5. Set the lancing device by extending it until it clicks into place. If it does not click, it is probably already in position from when you inserted the lancet.
6. The lancing device is ready to use. Do not lance your finger before the meter and test strips are ready to use.

**1) Consult the following for the suitable lancing depth: 1 - 2 for soft or thin skin, 3 - 4 for normal skin and 5 for thick or callous skin.**  
**2) Never use a lancet or a lancing device that belongs to another person. Sharing lancing devices and lancets may transmit blood borne pathogens, such as viral hepatitis.**

### Using the AST cap

It is generally recommended that a blood sample for a blood glucose test performed at home is taken from the fingertip. If you are not able to take a blood sample from your fingertip, you can also take one from another part of the body (AST), such as the ball of the hand or lower arm using the lancing device. In this case, exchange the protective cap on the lancing device with the AST cap. After inserting the lancet, put the transparent AST cap on the lancing device instead of the protective cap and tighten. Note that the AST cap is not intended to be used for a blood sample from the fingertip.

**We recommend, only to use the alternative site testing (AST), when:**

- in a pre-meal or fasting state (more than 2 hours since the last meal).
- at least 2 hours passed after the last meal
- at least 2 hours passed after the last intake of insulin and / or exercise

**Discuss the readings from an AST measurement with your doctor, if:**

- The glucose results are often fluctuating
- The results do not match the way you feel
- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.

**WARNING** In the case of low glucose (hypoglycemia), the blood sample must be drawn from the finger tip, as any changes in the blood glucose level are quicker to detect from a blood sample from the finger than from any other part of the body. The readings from a finger tip blood sample and from another part of the body can lead to the readings which radically differ from one another. Therefore, always consult your doctor before you perform a blood glucose test using a blood sample from another part of the body.

**AST Results:**

- If the blood glucose test result does not match how you feel, do a fingertip test to confirm the result again.
- Do NOT change your treatment just because of an alternative site result, do a fingertip test to confirm the result.
- If you often do not notice when your blood glucose is low, do a fingertip test.

**Caution:**

- Talk with your healthcare professional before you test on your forearm.
- Do NOT ignore symptoms of high or low blood glucose.
- Fingertip samples are able to show the rapid change of glucose faster than forearm samples.
- Do NOT change your treatment just because of a result.

The blood sample can also be taken from: (see illustration)


Proceed as follows:

1. Select a part of the body that is soft, not too densely covered with hair and not near a bone or a vein.
2. Massage the area gently to prepare the skin and to improve the circulation.
3. Hold the lancing device against the lancing spot for a few seconds and then press the trigger button.
4. Wait until a blood drop forms under the AST cap with a diameter of approx. 1.4 mm.
5. Remove the lancing device from the skin carefully and proceed in the same way as for the normal protective cap (see "Determining the blood glucose level").

**WARNING**

- Check with your healthcare professional before testing sites other than the fingertip.
- Do NOT use the first drop of blood sample when using the AST cap.

### Determining the blood glucose level

1. Insert a test strip into the meter in the direction of the arrow. The symbol appears automatically.
2. Press > to set (before a meal), (after a meal) or (before a meal with alarm after 2 hours) and press to confirm.
3. When the blood drop symbol flashes in the display, take a drop of blood from your finger tip. Massage the area gently to stimulate the blood circulation.
4. Place the lancing device on a finger tip (preferably at the side) and press the trigger button. Make sure, that the blood droplet does not smudge.


Apply blood to the edge of the test strip. Do not apply blood on top of the test strip.

Be sure to get enough blood on strip to make it to the confirmation window. Not enough blood on strip to get a test result.

5. Place the drop onto the blood sample area on the end of the test strip. Pay attention, that not any blood gets on the top of the test strip. To receive a correct test result, be sure to apply enough blood into the test strip's blood sample area - see the confirmation window as per drawing above.
6. When a sufficient amount of blood has been absorbed by the reaction cell you will hear a beep and in the display "OK" is shown.
7. The device will count down from 5 seconds (shown in display). After that, "OK" is shown and then the result appears in the display.

## Discarding used lancets

- Open the lancing device by twisting the protective cap in a clockwise direction and then remove it.
- Take out the used lancet (without touching it directly) and stick it into its protective cover.
- Move the sliding switch, located on the other side of the trigger button, which ejects the lancet. Dispose of the lancet in a container for biohazard material carefully to avoid harming other people.
- After discarding, wash your hands with soap and warm water. Rinse and dry thoroughly.

## Evaluating a test result

**WARNING**  
Never change the prescribed dose of medicine or treatment on your own initiative on the basis of one test result from a blood glucose test.

The **MediTouch® 2** blood glucose test strips work based on an improved technology (FAD-GDH) for a more exact and specific glucose measurement. They are calibrated to easily compare with laboratory test results.

Recommended premeal glucose concentration for nonpregnant adults with diabetes is between 80-130 mg/dL (4.44-7.22 mmol/L), and postmeal glucose concentration is < 180 mg/dL (9.99 mmol/L).

Consult your doctor about the range of blood glucose level valid for you.

## Unusual test results

If your test result is not what you expected, proceed as follows:

- Perform a control test, see section „Using the control solution“.
- Repeat the blood glucose test, see section „Determining the blood glucose level“.
- If your test result is still not what you think it should be, consult your doctor immediately.

**CAUTION**

- Extremely high humidity can affect the test results. A relative humidity of more than 90% can cause incorrect readings.
- If the proportion of red blood cells (haematocrit level) is very high (over 60%) or very low (under 20%), this could distort the test results.
- Studies have shown that electromagnetic fields can affect the test results. Do not perform a test near any devices which emit strong electromagnetic rays, such as microwaves, mobile phones etc.

\* Reference: *Diabetes Care* 2024;47(Suppl. 1):S111–S125 | <https://doi.org/10.2337/dc24-S006>

## Comparing your test results with a laboratory result

The question of how you can compare the blood glucose level of the meter with tests performed in the laboratory is frequently asked. Your blood glucose level can change quickly, especially after a meal, after taking medicine or strenuous activities.

Your blood glucose is affected by various factors and has different values at different times of the day. If you would like to compare the test result of your meter with a laboratory result, you must do the blood glucose test on an empty stomach. Therefore, it is advisable to do this in the morning. Take your meter with you to the doctor's and test yourself five minutes before or after a trained nurse has taken a blood sample from you. Take into account that the technology in the laboratory is different from your meter and that blood glucose meters for using at home generally produce slightly different results. To ensure the accuracy and precision of such important information, read the instructions included with the blood glucose test strips.

**NOTE**  
Make sure that you always record your test results with the date and time in your diabetes diary and label with the appropriate symbol for: before a meal ☉, after a meal ☉.

## Typical symptoms of high or low blood glucose

In order to better understand your test results, you can find some typical symptoms for high and low blood glucose herewith. In each case, you should contact your doctor about the therapy if you have noticed one of these symptoms.

### Result is higher than 240 mg/dL (13.33 mmol/L):

The test result is higher than reference normal range (80-130 mg/dL [4.44-7.22 mmol/L]).

### Possible symptoms may be:

Fatigue, increased appetite or thirst, frequent urination, blurred vision, headache, general aching, or vomiting.

### What to do:

- If you are experiencing any of these symptoms, test your blood glucose.
- If the result displayed is greater than 240 mg/dL (13.33 mmol/L) and you have symptoms of high blood glucose, contact your healthcare professional instantly.
- If the result does not match how you feel, follow the steps under "Unusual Test Results".

### Result is lower than 60 mg/dL (3.33 mmol/L):

The test result is lower than reference normal range (80-130 mg/dL [4.44-7.22 mmol/L]).

### Possible symptoms may be:

Sweating, trembling, blurred vision, rapid heartbeat, tingling, or numbness around mouth or fingertips.

### What to do:

- If you are experiencing any of these symptoms, test your blood glucose.
- If the result displayed is below 60 mg/dL (3.33 mmol/L) and you have symptoms of low blood glucose, contact your healthcare professional instantly.
- If the result does not match how you feel, follow the steps under "Unusual Test Results".

## Storing the test results

Your meter can save up to 480 test results, including the time and date. You can access the readings at any time. If the memory is full and you want to add a new test result, the oldest test result will be deleted automatically. Therefore, it is essential to enter the time and date correctly in your meter.

**NOTES**

- The content of the memory will not be deleted if you change the batteries. Just check that the time and date is still correct. You may have to reset the time and date after changing the batteries. To do this read „Setting the time and date“
- If the memory contains 480 test results and you want to add a new test result, the oldest test result saved will be deleted.
- Average values: Using the blood glucose readings from the previous 7, 14, 30 and 90 days, the MediTouch® 2 meter calculates the average value taken from the latest (480) to the first (001) test result and all those results AC (before a meal) and PC (after a meal) from the past 30 days.


## Accessing and deleting test results

You can access test results anytime without having to insert a test strip. Test results which are under certain criteria are marked with the appropriate symbol. When accessing saved test results, you can select according to these criteria by choosing the appropriate symbol:

before meals ☉

before meals with alarming 2 hours later ☉⚠

after meals ☉

- Press > to display the 30-day-average of the values saved under ☉. Further pressing of > displays the 30-day-average value of the saved values under ☉, then the total average values of the past 7/14/30/90 days. Press and hold ⏻ for 2 seconds to cancel this mode.
- Press ☉ to access the display mode for single values in memory. Press > to display the saved test results one after the other in sequence of 480 to 001 - in other words the latest entry will be displayed first and the oldest last.
- To delete a test result, simultaneously press and hold > and ☉ for 3 seconds. "dEL" appears in the display.
- Press ☉ to delete the result. The display shows „OK“.
- Press > to further display other saved values. Press and hold ⏻ for 2 seconds to cancel this mode. If no button is pressed, the meter switches off automatically after approx. 1,5 minutes.


If the following display appears during memory recall process, the device has no test results in memory so far. Perform a blood glucose test first, so that the device can save a test result.

## Integration in VitaDock® online

**medisana MediTouch® 2** features an integration option for transferring your measured values through the USB cable to the VitaDock app resp. online. Therefore you need the VitaDock software for your computer. Download the software from [www.medisana.com/software](http://www.medisana.com/software) as follows:

- Go to [www.medisana.com/software](http://www.medisana.com/software)
- Select "MediTouch 2" as your device.
- Download and install the VitaDock software on your computer.
- You will find a detailed instruction on the website, how to install and use the software.

- Connect the meter to your personal computer with the USB-cable. The - symbol appears in the display.
- The transfer of your measured values to the VitaDock® app resp. online will be started automatically if the software has been installed properly. Please check the page before for how to properly install it.
- After successful transfer, "OK" appears in the display and the device emits an acoustic signal.

## Maintenance and troubleshooting


### Display

Your meter needs little or no maintenance with normal use. It automatically tests its own systems every time you turn it on and lets you know if something is wrong. To ensure that the display is working properly, switch on the meter. Press and hold down the ☉-button for approx. 3 seconds so that you can see the whole display.

All the display elements must be clear to see with the accompanying symbol (please compare with the drawing on this page). If this is not the case, get in touch with the service centre.

**NOTES**

- A description of the symbols shown can be found at the beginning of this instruction manual.
- The device may also be switched on by inserting a test strip.


### Inserting Batteries

The meter uses two CR2032 Lithium Batteries. Batteries will normally last for more than 2000 tests. Other types of CR2032 Lithium Battery are also acceptable, but the capacity of test times may differ. Insert the batteries when you first use the meter or replace with new batteries when the „LP“ (low power) message and the low battery symbol appear on the display.

**1. The meter won't delete earlier records after replacing batteries.**  
**2. You should reset the time and date again after replacing the batteries.**  
**3. CR2032 Lithium Batteries are available at most stores.**  
**You may take the old batteries with you for replacement.**  
**4. Remove batteries when you will not be using the meter for one month or more.**

## Error displays

On no account should you change your medication due to a measurement, which may be incorrect. If in doubt, consult your service centre or your doctor.


A) What it means: Humidified / Used strips. The meter has detected a problem with the test strip.  
What to do: Repeat the test with a new strip.

B) What it means: Low batteries. The meter batteries do not have enough power to perform a test.  
What to do: Replace with new batteries.

C) What it means: Memory Error  
What to do: First replace the batteries. If error 005 appears again, get in touch with the service centre.


D) What it means: System Error  
What to do: First replace the batteries. If error 001 appears again, get in touch with the service centre.

E) What it means: Test result is higher than 630 mg/dL (35 mmol/L).  
What to do: Re-check your glucose level. If the result is HI again, obtain and follow instructions from your healthcare professional without delay.

F) What it means: The test result is lower than 20 mg/dL (1.1 mmol/L).  
What to do: This may require immediate treatment according to your healthcare professional's recommendations. Although this message could be due to a test error, it is safer to treat first and then do another test.

G) What it means: The test result is more than or equal to 251 mg/dL. The function of Ketone detector doesn't apply to the mode of control solution.  
What to do: Re-check your glucose level. If the result appears again, obtain and follow instructions from your healthcare professional without delay.

H) What it means: The "Hi" and thermometer icon appears: Temperature is too high, outside the required range of 10°C - 40°C (50°F - 104°F). This alerts users that an incorrect result may occur if the test continues.

The "Lo" and thermometer icon appears: Temperature is too low, outside the required range of 10°C - 40°C (50°F - 104°F). This alerts users that an incorrect result may occur if the test continues.

What to do: Relocate the meter to a location with temperature between 10°C - 40°C (50°F - 104°F).

I) What it means: Volume detector error. The volume of blood or control solution is NOT enough.  
What to do: Replace the batteries first. If **Err** appears again, please contact your local dealer.

J) What it means: No result in memory. The meter doesn't run any test yet.  
What to do: You can still perform a blood glucose test and get an accurate test result.

## Cleaning and maintenance

### Cleaning Your Meter

Caring for your MediTouch® 2 SMBG system does not require special cleaning. Please keep the meter free of dirt, dust, bloodstain and water stains.

Follow these guidelines carefully to help you get the best performance possible:

- Make sure the meter is turned off.
- Gently wipe the meter's surface with a soft cloth slightly dampened with ethanol (70–75%).
- DO NOT let any moisture get into the test strip slot.
- DO NOT spray any cleaning solution directly onto the meter.
- DO NOT put the meter under water or liquid.
- DO NOT pour liquid into the meter.


### Cleaning Your Lancing Device

- To clean the lancing device, wipe it with a soft cloth dampened with water and mild detergent.
- DO NOT place the entire device under water.
- To disinfect the cap after cleaning, place it in 70%-75% rubbing alcohol for 10 minutes at least once a week. Allow the cap to air-dry after disinfecting.

## Reset the meter

To reset the meter (Attention: All saved data will be deleted!), open the cover of the battery compartment on the back side of the unit.

Press and hold the reset button for approx. 3 seconds to reset the device.

## Technical specifications

Name and model:	<b>medisana</b> Blood glucose meter <b>MediTouch® 2 / Model: MDT2</b>
Measuring method:	electrochemical biosensor technology
Measuring range:	20 - 630 mg/dL (1.1 - 35.0 mmol/L)
Test time:	approx. 5 seconds
Memory:	480 test results with time and date
Operating conditions:	Temperature 10°C – 40°C (50°F – 104°F), relative humidity 40 - 85 %
Storage / transport conditions:	Temperature 2°C – 30°C (35.6°F - 86°F), relative humidity 40 - 85 %
Blood sample volume:	0,6 µL
Sample material:	Fresh blood from finger tip, palm of the hand or arm (capillary whole blood)
Haematocrit value (Htc):	20 – 60 %
Power supply:	2 x 3V CR2032 lithium batteries
Battery operation:	over 2,000 tests
Automatic switch-off:	after approx. 1.5 minutes
Dimensions (display):	1.46" x 2.20" (37.0 x 56.0 mm)
Dimensions (base unit):	1.97" x 3.86" x 0.39" (50.0 x 98.0 x 10.0 mm)
Weight:	1.06 oz (30.0 grams) w/o batteries
Article number:	79030
EAN number:	40 15588 79030 0

**CE 0537**

**In accordance with our policy of continual product improvement, we reserve the right to make technical and visual changes without notice.**

Electromagnetic compatibility: The device complies with the EN 60601-1-2 standard for electromagnetic compatibility. Enquire at medisana for details on this measurement data.

The meter is certified in accordance with the requirements of Directive 98/79/EC for in vitro diagnostic devices.

**Accessories**

- Enquire at your local supplier or service centre:
- 100 MediTouch® 2 lancets Art.-Nr. 79028
- 1 MediTouch® 2 lancing device Art.-Nr. 79002
- 1 MediTouch® 2 control solution Art.-Nr. 79038
- 50 MediTouch® 2 test strips Art.-Nr. 79038

**Disposal instructions**

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. **Teststrips and lancets:** Always dispose of the test strips and lancets in a way that prevents injury and the spread of infection to others. Consult your local authority or your supplier for information about disposal.

**Comparison of Measured Values With Laboratory Value - Precision**

Three lots of the MediTouch 2 blood glucose test strips have been tested to assess the precision of blood glucose measuring system. This includes a repeat assessment using venous blood and a laboratory precision assessment using the control material. The blood glucose content of the venous blood samples ranges from 42.7 to 418.0 mg/dL and control material from three concentrations is used. Results of the repeat precision measurements:

Sample	Venouse blood mg/dL (mmol/L)	Grand mean value mg/dL (mmol/L)	Pooled standard deviation mg/dL (mmol/L)	Pooled SD with 95% CI mg/dL (mmol/L)	Pooled coefficient of variation (%)
1	43 (2.4)	36 (2.0)	2.0 (0.1120)	2.01~2.02 (0.1118~0.1122)	5.6
2	62 (3.4)	59 (3.3)	3.5 (0.1933)	3.47~3.49 (0.1928~0.1939)	5.9
3	121 (6.7)	127 (7.1)	4.1 (0.2282)	4.06~4.16 (0.2256~0.2237)	3.2
4	201 (11.2)	214 (11.9)	6.7 (0.3712)	6.59~6.78 (0.3659~0.3764)	3.1
5	317 (17.6)	330 (18.3)	10.1 (0.5607)	10.05~10.15 (0.5580~0.5634)	3.1
6	418 (23.2)	433 (24.1)	14.5 (0.8030)	14.36~14.58 (0.7970~0.8091)	3.3

Results of the intermediate precision measurement:

Sample	Grand mean value mg/dL (mmol/L)	Pooled standard deviation mg/dL (mmol/L)	Pooled SD with 95% CI mg/dL (mmol/L)	Pooled coefficient of variation (%)
1	71 (3.9)	1.012 (0.0562)	1.010~1.014 (0.0561~0.0563)	1.4
2	136 (7.6)	1.431 (0.0794)	1.418~1.443 (0.0787~0.0801)	1.1
3	351 (19.5)	2.823 (0.1567)	2.819~2.826 (0.1565~0.1569)	0.8

**System Accuracy**

The MediTouch 2 blood glucose meter in comparison with the YSI. Three lots of MediTouch 2 Blood Glucose Test Strips have been tested to assess the system accuracy of the MediTouch 2 blood glucose measuring system and to compare it with the reference method in which capillary whole blood concentrations of 34.4 to 442.8 mg/dL have been used.

System accuracy results for glucose concentrations <100 mg/dL (<5.55 mmol/L):

within±5mg/dL (Within ± 0.28 mmol/L)	within±10mg/dL (Within ± 0.56 mmol/L)	within±15mg/dL (Within ± 0.83 mmol/L)
61/186 (32.8 %)	117 /186 (62.9 %)	181/186 (98.4 %)

System accuracy results for glucose concentrations ≥ 100 mg/dL (≥5.55 mmol/L):

within ±5%	within ±10%	within ±15%
205/414 (49.5 %)	339/414 (81.9%)	398/414 (96.1%)

System accuracy results for glucose concentrations between 34.4 mg/dL (1.91 mmol/L) and 442.8 mg/dL (24.58 mmol/L):

Within ± 15 mg/dL or ±15% (Within ± 0.83 mmol/L or ±15%)
579 / 600 (96.5 %)

In comparison to the YSI, the MediTouch 2 met the EN ISO 15197:2015 standard, whereby 95% of the blood glucose values measured have to fall within the following zones: either ±0.83 mmol/L (±15 mg/dL) of the measured average value when using the reference measuring procedure for blood glucose concentrations <100 mg/dL (<5.55 mmol/L) or ±15% for blood glucose concentrations of ≥ 100 mg/dL (≥5.55 mmol/L). 99% of the individual measured blood glucose values must fall within zones A and B of the Consensus Error Grid (CEG) for diabetes type 1.

**Performance evaluation by the user**

A study to assess the glucose values of blood samples of capillary blood from the fingertips, which were obtained from 106 individuals that had no special training, produced the following results: 96.9% within ±15mg/dL (± 0.83 mmol/L) and 96.2% within ± 15% of the values obtained in the medical laboratory with glucose concentrations of at least 100 mg/dL (5.55 mmol/L). You will find further details and information regarding blood glucose results and various technologies in generally relevant specialist medical literature.

**Limitation**

The test strips are used for fresh capillary whole blood samples.

1. DO NOT use serum or plasma sample.
2. DO NOT use anticoagulant NaF or potassium oxalate for venous sample preparation.
3. DO NOT use neonate blood sample.
4. Extreme humidity may affect the results. A relative humidity greater than 90 % may cause incorrect results.
5. The system should be used at temperatures between 10°C and 40°C (50°F and 104°F). Outside this range, the system may get incorrect results.
6. DO NOT reuse the test strips. The test strips are for single use only.
7. Hematocrit: The hematocrit between 20% and 60% will not affect the results. Hematocrit below 20% may cause higher results. Hematocrit above 60% may cause lower results.
8. Altitude up to 3,048 meters above sea level has no effect on readings.

**Healthcare Professionals – Please note these additional Limitations**

1. If the patient has the following conditions, the result may fail:
  - ◆ Severe dehydration
  - ◆ Severe hypotension (low blood pressure)
  - ◆ Shock
  - ◆ A state of hyperglycemic-hyerosmolar state (with or without ketosis)
2. Lipemic samples: Cholesterol level up to 400 mg/dL (10.32 mmol/L) and triglycerides up to 800 mg/dL (9.04 mmol/L) do not affect the results. Grossly lipemic patient samples have not been tested and are not recommended for testing with MediTouch 2 Glucose Meter.
3. Critically ill patients should not be tested with MediTouch 2 glucose meters.
4. DO NOT use during xylose absorption testing. Xylose in the blood will interfere Self-Monitoring Blood Glucose System.
5. Interfering Substances depend on the concentration. The below substances up to the test concentration will not affect the test results.

Concentrations of the interference tested	Bias   Glucose level		
	50-100 mg/dL (2.8-5.6 mmol/L)	250-350 mg/dL (13.9-19.4 mmol/L)	
Acetaminophen	7 mg/dL (0.46 mmol/L)	8.1 (0.45)	5.3 %
Ascorbic Acid	4 mg/dL (0.26 mmol/L)	6.6 (0.37)	5.8 %
Bilirubin-unconjugated	3.3 mg/dL (0.056 mmol/L)	0.2 (0.01)	5.2 %
Cholesterol	400 mg/dL (10.32 mmol/L)	9.6 (0.53)	7.2 %
Creatinine	30 mg/dL (2.7 mmol/L)	1.3 (0.07)	1.6 %
Dopamine	2.2 mg/dL (0.14 mmol/L)	8.0 (0.44)	3.2 %
Galactose	20 mg/dL (1.11 mmol/L)	6.2 (0.34)	2.5 %
Gentisic Acid	7 mg/dL (0.45 mmol/L)	9.8 (0.54)	3.6 %
Glutathione	1 mg/dL (0.03 mmol/L)	3.7 (0.21)	6.5 %
Haemoglobin	300 mg/dL (0.05 mmol/L)	3.8 (0.21)	5.2 %
Ibuprofen	50 mg/dL (2.43 mmol/L)	3.9 (0.22)	2.7 %
Icodextrin	1094 mg/dL (0.66 mmol/L)	5.4 (0.30)	4.8 %
L-Dopa	2 mg/dL (0.10 mmol/L)	10.0 (0.56)	8.7 %
Maltose	278 mg/dL (7.78 mmol/L)	2.7 (0.15)	4.4 %
Methyldopa	4 mg/dL (0.19 mmol/L)	9.0 (0.50)	3.7 %
Pralidoxime Iodide	5 mg/dL (0.14 mmol/L)	2.8 (0.16)	3.3 %
Sodium Salicylate	40 mg/dL (2.5 mmol/L)	4.3 (0.24)	2.2 %
Tolbutamide	100 mg/dL (3.70 mmol/L)	1.4 (0.08)	2.3 %
Tolazamide	2.5 mg/dL (0.08 mmol/L)	2.5 (0.14)	3.6 %
Triglycerides	800 mg/dL (9.04 mmol/L)	9.3 (0.52)	5.6 %
Uric acid	16.5 mg/dL (0.99 mmol/L)	7.2 (0.40)	4.0 %
Xylose	9.5 mg/dL (0.63 mmol/L)	7.0 (0.39)	7.5 %

**Lancets**

Correspond to the REGULATION (EU) 2017/745

Manufacturer:  
**Shandong Lianfa Medical Plastic Products Co.,Ltd.**  
 No.1 Shuangshan Sanjian Road, 250200 Zhangqiu City, Jinan, Shandong, P.R. China  
 Tel: +86 531 6132-8777 linkfar@sdliafa.com

**EC REP**

**Linkfar Healthcare GmbH**  
 Nierrehinstraße 71, 40474 Düsseldorf, Germany  
 TEL.: +49-21130036618

**Lancing device**

**GMMC**  
 Room 1112, Ace Tower 9th Bldg., 130, Digital-ro, Geumcheon-gu, Seoul, South Korea  
 Tel : +82 2 2104 0470  
 info@gmmcgroup.com

**EC REP**

**GMMC, S.L**  
 Calle Carolina Alvarez, 41-1, Valencia, Spain  
 Tel: +34 636 170 691

**Warranty and repair terms**

Your statutory warranty rights are not restricted by our guarantee below. Please contact your supplier or the service centre in case of a claim under the warranty. If you have to return the device, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **medisana** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Where defects in materials or workmanship arise we will repair or replace free of charge with in the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - b. Spare parts which are subject to normal wear and tear, such as battery, lancing device and disposable items etc.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

**EPS BIO Technology Corp.**  
 No.8, R&D RD. III, Hsinchu Science Park, Hsinchu, Taiwan 30076

**EC REP**

Obelis s.a.  
 Boulevard Général Wahis 53, B-1030 Brussels, Belgium



**medisana GmbH**  
 Carl-Schurz-Str. 2, 41460 NEUSS, Germany

**Service**

You can find out which is your responsible service centre from your supplier.

The current version of this instruction manual can be found under [www.medisana.com](http://www.medisana.com)

- 19 Time (hours:minutes)  
AM = before noon; PM = after noon
- 20 Memory slot
- 21 Blood glucose reading
- 22 Unit of measure
- 23 Battery symbol (low battery)
- 24 Ambient temperatur error
- 25 Insert test strip
- 26 Apply blood or control solution
- 27 Active USB connection
- 28 System error
- 29 Warning regarding a possible diabetical ketoacidosis.  
Consult doctor!
- 30 Control solution measurement mode
- 31 Date (month:day)
- 32 : Alarm (before meals with automatic alarm after 2 hours to remind you to perform a test after meals)
- Before meals (AC)
- After meals (PC)

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