DE/GB/NL/FR/ES/IT/PT/GR/FI/SE/NO/DK

medisana®

OL 800 /
OL 850

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GB Heated Pad

NL Warmtedeken

FR Couverture chauffante

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IT Coperta termica

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GB Safety instructions

Read the instructions for use carefully, especially the safety instructions, before using the appliance and keep the instructions for future use. These brief instructions are an integral part of the instructions for use. The complete instructions for use are available for download at <https://www.medisana.com>. If you pass the unit on to a third party, be sure to include the instructions for use.

Power supply: Before connecting the unit to the power supply, make sure that it is switched off and that the specified electrical voltage is observed. • Always switch off the appliance after use with the button on the control unit and disconnect the power supply. • Never carry, pull or turn the device by the cable. • The cable and control unit of the appliance can give rise to risks of entanglement, strangulation, tripping or treading if not correctly arranged. The user shall make sure that excess ties and electric cords shall be arranged in a safe way. • *For special people:* Do not use the device on children or persons who are disabled, asleep or insensitive to heat (persons who cannot react on overheating). • This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. • Children shall not play with the appliance. • Children under the age of three are not allowed to use this appliance due to their inability to react to overheating. • The appliance is not allowed to be used by young children over the age of three unless the controls have been pre-set by a parent or guardian, or unless the child has been adequately instructed on how to operate the controls safely. • If you have any health concerns, talk to your doctor before use. • *Operating the device:* Only use the device for its intended use as in the instruction manual. • If misused, all guarantee rights shall become null and void. • Before each use, carefully check the article for signs of wear and/or damage. Do not operate the appliance if you notice wear, damage or signs of improper use. • If an error occurs, check that the power supply is properly established. Do not open the unit. There are no user-serviceable parts inside the unit. • The heated blanket must not be operated folded or bent. • Do not immerse the device in water or other liquids. If liquid does enter the device, disconnect the power supply immediately. • Do not use the device in a wet environment (e.g. in the bathroom or similar). • Do not keep or place any electrical device on a surface where it is in danger of falling into a bathtub or sink. • Do not use the device while you are sleeping or lying in bed. Danger of burns or heat stroke! • The device is not suitable for medical use (e.g. in hospitals). • Very long periods of use can cause overheating. Always allow the appliance to cool down before using it again. • Never operate the appliance unattended. • Do not attach safety pins or other pointed or sharp objects to the appliance or prick it. • Do not cover the unit when it is on. Never use it under blankets or pillows. There is a risk of fire, electric shock and injury. • Never store or use the appliance near an electric cooker or other heat sources. • *Maintenance and cleaning:* Do not repair the unit yourself in the event of a fault. Repairs may only be carried out by an authorised dealer or a suitably qualified person. • Cleaning and user maintenance shall not be made by children without supervision. • When storing the appliance, allow it to cool down before storing it. • Do not place any objects on top of the item during storage to avoid sharp kinks. • Do not crease the appliance by placing items on top of it during storage. • Examine the appliance frequently for signs of wear or damage. If there are such signs, if the appliance has been misused or does not work, contact the customer service before switching it on again.

❶ On/Off button; ❷ USB plug for power bank*; ❸ Integrated pocket for powerbank*; ❹ Strap to roll up

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