medisana
 Carl-Schurz-Str. 2,
 41460 NEUSS, Germany

REF 61134
EAN 40 15588 61134 6

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Heating pad HP 519 GB	
1 Control unit HD20-A	6 Plug connection
2 On/Off button	7 Button to change heating level
3 Cable with plug	8 LED display of the heating level (available levels: 1 ... 4)
4 Heating pad	

Safety Instructions
Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. These brief instructions are an integral part of the instructions for use. The complete instructions for use are available for download at https://www.medisana.com. Should you give this device to another person, it is vital that you also pass on the instructions for use.

- Please ensure that the packing is kept away from the reach of children! Risk of suffocation!
- The device must only be used for its intended purpose as described in the instruction manual. Using the appliance for any other purpose invalidates the warranty.

- The device is not to be used for commercial use or medical purposes.
- When the device is likely to be used for an extended period, for instance by the user falling asleep, controls should be adjusted to a setting recommended for continuous use.
- Check the article carefully for signs of wear and/or damage before each use. Do not operate the device if you notice wear, damage or signs of improper use, or if it does not function correctly. Instead, contact an authorised service centre in such a case.
- Before connecting the heating pad to the power supply, please check that it is switched off and that the mains voltage of your power supply corresponds to the voltage listed on the label.
- Keep mains cable and control unit away from heat, humidity and liquids. Never touch the mains plug when your hands are wet or when you are standing in water.
- Only operate the heating pad with the corresponding switching unit HD20-A.
- Do not attempt to touch a heating pad which has fallen into water. Pull the plug out of the socket immediately.
- To disconnect the heating pad, always pull the mains plug out of the mains socket. Never pull on the cable!
- Do not carry, pull or twist the heating pad with the mains cable.
- Make sure that the power socket and plug are easily accessible so you can unplug the mains quickly if necessary.
- The cable and control unit of the appliance can give rise to risks of entanglement, strangulation, tripping or treading if not correctly arranged. The user must ensure (e.g. using cable ties) that all electrical cables are arranged in a safe way.
- The device complies with the guidelines on electromagnetic compatibility. If you wear a pacemaker, this could be affected. Consult your doctor and the manufacturer of your pacemaker before using this device.
- If you have any medical concerns, please discuss using the heating pad with your doctor before using it.
- Do not use the device on children or persons who are disabled, asleep or insensitive to heat (persons who cannot react on overheating).
- Children below the age of 3 years are not allowed to use this device, because they are not able to react on overheating.
- The device must not be used by young children over the age of 3 years, unless the control unit has been adjusted accordingly by its parents or legal guardians and unless the child has been well instructed how to use the control unit properly.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- The appliance can be used by children older than 3 years and younger than 8 years under supervision and with the control always set to minimum temperature value.
- Children shall not play with the appliance.
- Please ensure that the polythene packing is kept away from the reach of children! Risk of suffocation!
- Do not use the heating pad on parts of the body which are swollen, inflamed or injured.
- Consult your doctor, if you have pain in muscles or joints for a prolonged period of time.
- Only use the heating pad as intended and in accordance with the instruction manual.
- The warranty will become void if it is used improperly.
- Do not use the unit for animals.
- This heating pad is intended only for private use and not for commercial purposes or for use in hospitals.
- The heating pad is subject to aging. Check heating pad, mains cable and control unit carefully for damage before each use. A faulty unit must not be used.
- Do not use the appliance if it is damaged, if its power cord is damaged, if it is not in perfect working order, or if it has been dropped or fallen into water. Send the unit to the service centre for repair to avoid potential hazards.
- Keep the unit out of rain and do not let it get wet. It must not be used out-of-doors.
- Never place the unit immediately next to a heater or radiator.
- Ensure that the unit does not come into contact with pointed or sharp objects and do not push pins into it.
- The heating pad must not be used if it has been folded, bent or crushed.
- Do not use the heating pad underneath quilts, blanket or pillows or cushions.
- Do not use the heating pad when wet and only use it in dry environments (not in the bathroom or similar environments).
- The control unit must not be placed on or under the heating pad during operation and it must not be covered.
- Do not sit on the heating pad, but place it like a cushion on the body part, as required.
- Do not leave the heating pad unattended. Always switch off the heating pad when not in use and pull the plug out of the socket. Never use a time switch.
- Do not fall asleep while the heating pad is switched on and the plug is inserted into the socket.
- Ensure that you follow the safety notes, as prolonged use at a high setting can result in burns to the skin. Check the skin areas under the heating pad regularly to prevent burns or blisters.
- If use of the heating pad is unpleasant or painful, discontinue using it immediately.
- You may only clean the heating pad. To avoid potential hazards, you must never attempt to repair it yourself. Contact the service centre.
- Cleaning and user maintenance shall not be made by children without supervision.
- Before cleaning the unit, switch it off and unplug it from the mains. Never immerse the control unit in water or other fluids.
- Let the heating pad cool down completely before storing it. Store it in a place that is dry and free from dust, preferably in the original packaging, until you use it next.
- Do not create the appliance by placing items on top of it during storage.
- If a fault occurs, do not attempt to repair the unit yourself. Repairs must only be carried out by an authorized specialist dealer or other suitably qualified personnel.
- A damaged power supply cable must be replaced with a power cable of the same type.

Using the device
 Make sure, the connection cable of the control unit is properly connected to the plug of the heating pad. Plug the unit into the mains outlet and press the On/Off button ➊. The device starts at lowest level (level 1), which is displayed at ➋.

After a few minutes you will be able to feel the pad getting noticeably warmer. If you require a higher temperature setting, press the Heat button ➌ several times until you reach the required level. (The maximum level is level 4, after that the level returns to the lowest level 1.) If the heating pad gets too hot for you, press the Heat button ➍ several times until you reach the required level. After a short time, the heating pad will cool down. The heating pad may cool down after a short period of time. This is normal. The fast re-heating according to the setting will follow immediately.

For continuous use, you have to always select the lowest level 1.

To switch the heating pad off, press the On/Off button ➋. The LEDs ➌ will go out, indicating that the unit is switched off. Unplug the unit from the mains outlet.

The heating pad will switch off automatically after approx. 90 minutes of continuous use (unless you have manually switched it off before). Unplug the unit from the mains outlet once you have finished using it.

- Cleaning and maintenance**
- Before cleaning the pad, unplug it from the mains outlet and let it cool down for at least ten minutes.
 - The heating pad is equipped with a removable control unit (HD20-A). Unplug the control unit from the heating pad by pulling out the connection cable from the plug on the heating pad.
 - You may wash the heating pad gently by washing machine at 30°C.
 - However, the best way is to put the heating pad in a bath with lukewarm water and some mild detergent and then gently squeeze it out. Rinse the heating pad several times to remove all the detergent.
 - Never use strong detergents or hard brushes.
 - To dry, lay the pad flat on a moisture-absorbent towel or mat.
 - Only use the blanket again once it is completely dry.
 - Reconnect the control unit to the heating pad in the correct manner.
 - You should untangle the cable from time to time.
 - Store the unit in a clean and dry place, laid out flat with nothing covering or resting on top of it.

Voltage supply: 220 – 240V~ 50/60Hz approx. 80W

Disposal: This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

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