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Sensor- spitze 1 ON/OFF-Taste (Start) 3

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**GB** LCD display 2 Cap 4

Sensor tip 1 ON/OFF-button (Start) 3

**ORALE MESSUNG:** Place the thermometer with the sensor tip 1 under the tongue and close the mouth. Breathe through the nose during measurement to prevent the reading from being distorted. **IN THE RECTUM:** Many patients find rectal measurement unpleasant, but it does deliver accurate results. Carefully insert the sensor tip 1 approx. 3.5 cm into the rectum. **MEASUREMENT UNDER THE ARM:** Clean the armpit before taking the temperature. Place the thermometer with sensor tip 1 in the armpit. Lower the arm and lightly press it against the body during measurement.

**SAFETY INSTRUCTIONS**  
**Important notes! Be sure to keep!**  
 Please read these instructions for use carefully before using the device, in particular the safety instructions and keep the instructions in a safe place for future reference. If you give the device to someone else, be sure to include these instructions for use.  
 The thermometer contains small parts that could be swallowed by children. Do not leave the thermometer unattended in the hands of children.  
 Use the thermometer only for taking body temperature under the tongue, in the armpits, or in the rectum. Do not use the thermometer in the ears.  
 Use the thermometer only for measuring body temperature. The thermometer is not suitable for professional use. Do not attempt to repair the thermometer yourself. If it no longer functions properly, do not continue to use it. Ask your dealer and have repairs done only by **medisana** service.  
 Observe the operating and storage conditions in order to keep the **medisana** FTC thermometer fully functional for a long time.  
**SCOPE OF DELIVERY**  
 Please first check whether the device is complete. The contents include:  
 - the **medisana** FTC thermometer  
 - Battery type AG3 or LR41 (already installed)  
 - Storage box  
 - Instructions for use

**Remove/insert battery**  
 The **medisana** FTC thermometer is ready to use out of the package as the battery is already pre-inserted. However, if after a certain period of operation, a battery symbol appears on the display, the battery must be replaced. To do this, pull off the cap 4 of the battery compartment and remove the battery compartment about 18 mm from the thermometer. Use a suitable tool such as a ballpoint pen tip to loosen the button cell battery carefully. Insert the new button cell battery with the "+" symbol pointing upwards. Slide the battery compartment back into the thermometer. Close the battery compartment again with the cap.  
 Range of life: 32,0°-42,0°C  
 Battery of measurement: Measurement accuracy: ± 0.1°C (deviation from the actual body temperature in the measuring range of 35,5°-42,0°C) not extrapolating  
 Storage and transport conditions: 10,0°-40,0°C  
 Dimensions (L x W x H): ca 10 mm x 128 mm x 19 mm  
 Weight: ca 9.8 g  
 EAN Code: 40 15588 77030 2

**Faking temperatures: fast and error-free!**  
 Fevers are often associated with fatigue, body aches, shivering, and chills. The course of the disease can best be controlled by measuring the fever two to three times a day at intervals of 4 to 6 hours. The measurement should always be performed in the same area of the body and at about the same time for the results to be comparable. The **medisana** FTC thermometer makes it quick and easy to take a current body temperature. Frequent errors such as not leaving the thermometer in place long enough can almost be ruled out. A temperature measurement can be carried out in the following areas: a) under the tongue (oral), b) in the armpit (axillary), c) in the rectum (rectal).

**Measurement**  
 Turn the thermometer on by pressing the Start button 3. All of the symbols will appear for about two seconds, followed by a short beep. The most recent reading will then be briefly displayed on the LCD display 2 (with a small M appearing on the right edge of the display). If the temperature is less than 32°C, you will see "Lo" and a flashing "1" on the LCD display 2. This means the **medisana** FTC thermometer is ready to use. Position the sensor tip under the tongue, in the armpit, or in the rectum. The thermometer will then automatically start measuring the temperature. During the measurement, the current temperature is displayed continuously. As soon as the measurement is completed, 10 beeps sound and the measured temperature is shown on the display. If the thermometer detects a temperature of more than 37,7 °C, a fever alarm will sound in the form of 10 triple beeps. The thermometer switches off automatically eight minutes after the measurement. To extend battery life, however, you should turn the thermometer off by pressing the Start button 3 once you have read the result.

**Taking the temperature in different parts of the body**  
 The **medisana** FTC thermometer is suitable for taking temperatures under the tongue, in the armpit, or in the rectum. It is completely waterproof. After about 1 minute, 10 short beeps will sound and the measurement result is shown on the LCD display 2. Then the thermometer can be removed. To ensure a safe and accurate measurement, you should consider the following when positioning the thermometer:

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