medisana®

FS 881

Manual

Service

medisana
Carl-Schurz-Str. 2,
41460 NEUSS, Germany

REF

88380

EAN

40 15588 88380 4

Fußsprudelbad FS 881

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Use

Fill the device with warm water. Place the device on a firm floor area and then plug in the power supply to the mains. Sit down comfortably in front of the device and switch on the device by pressing the function button ①. Your feet will be stimulated by vibration of the massage nodes ② at the bottom of the device and bubbles. You may use the acupuncture point massage area ③ to freely apply to desired parts of your feet, too.

Rolling your soles over this massage area increases blood circulation and has an invigorating effect on the feet.

After finishing the footbath, press the function button ① again to switch the device off. Pull the plug out of the mains supply and empty the foot spa. To empty, tip the unit, allowing water to pour out along the back drain spout. Make sure, no water will come in contact with the function button ①.

Power supply: 220-240V~50Hz

Disposal: This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your local authority or your supplier for information about disposal.

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• We advise consultation with your doctor if you are pregnant.

• Stop a massage if you experience pain during the application or swellings crop up after use.

• The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.

• Only use the device for its intended use as described in the instruction manual. If misused, all guarantee rights shall become null and void.

• The device is not determined for commercial use, but only for domestic use.

• Do not use the device outdoors or in damp rooms.

• Make sure that you use your foot spa on a safe and even surface that is not sensitive to moisture.

• Only use the device with water and not with other liquids.

• Only use appropriate bath additives – do not use foaming bath additives or bath salts.

• Do not stand in the device as it is not designed to take the burden of full body weight.

• Never cover the device during operation and do not use it under blankets or cushions.

• If water leaks from the appliance, the appliance should no longer be used.

• Do not leave the Foot Spa unattended while it is operating.

• Always check the temperature of the water before placing your feet into the Foot Spa.

• The Foot Spa cannot be used to warm cold water. You must fill the unit with warm water properly as you need.

• Never fall asleep while using the Foot Spa.

• Never drop or insert any object into any opening.

• Do not operate where aerosol (spray) products are being used or where oxygen is being administered.

• In case of a malfunction, do not try to repair the device yourself. Not only does the warranty expire, but serious danger can result (fire, electric shock, injuries). Only authorised service offices should carry out a repair.

• Cleaning and user maintenance shall not be made by children without supervision.

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