medisana®

AC 850

- DE Cellulite-Massagegerät
- GB Cellulite Massager
- NL Cellulite massageapparaat
- FR Appareil de massage anti-cellulite
- ES Masajeador para celulitis
- IT Massaggiatore anticellulite
- PT Aparelho de massagem anticelulite
- PL Masażer antycellulitowy
- CZ Masážní přístroj proti celulitidě
- HU Cellulitisz masszáz-szélkülék
- GR Συσκευή μασάζ κατά της κυτταρίτιδας
- TR Selülit Masaj Cihazı
- FI Seluliitin hierontalaite
- SE Massageapparat mot celluliter
- RU Антицеллюлитный массажер

PDF online

www.medisana.com

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Safety Instructions

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. These brief instructions are an integral part of the instructions for use. The complete instructions for use are available for download at <https://www.medisana.com>. Should you give this device to another person, it is vital that you also pass on the instructions for use.

- The device is intended for indoor private use only and should not be employed for commercial purposes.
- Only use the device for the according to the instructions for use. The guarantee is invalidated if the device is misused.
- Avoid contact with sharp or pointed objects.
- Do not use this device if you are pregnant.
- Persons with a pacemaker, artificial joints or electronic implants should always consult a doctor before using the device.

- You must not use the device if you suffer from one or more of the following diseases or medical conditions: Circulatory problems, varicose veins, open wounds, diabetes, fever, heart problems, bruises, skin tears or other skin injuries, phlebitis, tumours, thrombosis, osteoporosis, spinal fractures or slipped discs.
- Do not use the appliance on the following parts of the body: Head, face, front of neck, chest, abdomen, bones or joints (e.g. elbows, knees, etc.), kidney and genital area.
- Excessive pressure when massaging with the device can lead to injuries!
- Special care is necessary when the device is used on, by or near children, sick and helpless persons.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- If you feel any pain from use of this device, or find the massage unpleasant, stop the treatment and consult your doctor. Keep the appliance away from heat, hot surfaces, moisture and liquids. Never touch the appliance with wet hands or when standing in water. Do not reach for a unit that has fallen into water.
- Do not cover the device when it is switched on. Never use it under blankets or pillows. There is a risk of fire, electric shock and injury.
- Never store or use the appliance near an electric stove or other heat sources.
- Do not sit or stand on the unit.
- Do not use the appliance together with other appliances.
- Do not use the appliance before going to bed. The massage has a stimulating effect and can delay falling asleep.
- Do not fall asleep while using the massager on yourself. Using it on the same spot for too long can cause discomfort and burns.
- Do not repair the device yourself in case of a malfunction. Repairs may only be carried out by an authorized dealer or an appropriately qualified person.
- Check the unit for damage before use. A defective unit must not be put into operation!
- Cleaning and user maintenance shall not be made by children without supervision.
- Switch off the device after each use.
- Do not use the device immediately after eating or drinking alcohol or when you are tired.
- Do not drop the appliance or subject it to strong shocks.
- Do not apply the appliance while operating vehicles or machinery.
- Make sure that the cable does not become a tripping hazard. It must not be kinked, pinched or twisted. If the cable is defective, it must not be used any longer.
- A damaged mains cable can be replaced with a mains cable of the same type.
- If you have long hair, it is advisable to fix the hair appropriately behind or on top of the head so that there is no risk of tangling when using the massager.

Cleaning and care

- For cleaning, make sure that the adapter with cable is not connected and let the unit cool down for at least 10 minutes after use.
- Never use aggressive cleaning agents, thinners, petrol, strong brushes or similar. The surface of the unit is sensitive.
- The unit should be cleaned with a soft, slightly damp cloth, possibly with a mild soap solution. To make cleaning easier, you can remove the massage rollers and clean them separately. Do not mount the rollers again until they are completely dry.
- Clean the adaptor with cable with a clean, slightly moistened cloth.
- Never immerse the device or the adaptor in liquids.
- Only use the unit again when it is completely dry and make sure that no moisture can get into the unit!
- It is best to store the unit in its original packaging and keep it in a

medisana GmbH
Carl-Schurz-Str. 2, 41460 NEUSS,
Germany

REF 88540 EAN 40 15588 88540 2

clean and dry place out of direct sunlight.

Power supply: Adapter with connection cable, Input: 100–240V~ 50/60Hz, 400mA | Output: 12V=, 1200mA

Disposal: This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

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