**PDF online**  
 www.medisana.com

medisana GmbH  
 Carl-Schurz-Str. 2, 41460 NEUSS,  
 Germany

REF 88963 EAN 4015588 88963 9

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- 1 On/Off button
- 2 Button to select between 2 massage directions
- 3 Button to select between 3 massage intensity levels
- 4 HEAT button, switches heat function on/off
- 5 Hand strap (each 1 on right and 1 on left side)
- 6 Light indicator (LED)
- 7 Removable neck cover

### Safety Instructions

**Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. These brief instructions are an integral part of the instructions for use. The complete instructions for use are available for download at <https://www.medisana.com>. Should you give this device to another person, it is vital that you also pass on the instructions for use.**

### About the power supply

- Only charge the unit's battery using the USB charging cable supplied. Never leave the unit unattended when it is connected to the mains (for charging) or in operation.
- Never charge the device in a flammable environment (e.g. flammable objects, gasoline, flammable gases, paint, etc.)
- Only connect the unit for charging so that the cable and power connector are freely accessible.
- The built-in batteries are not intended for replacement or removal by the end user. If problems occur with the power supply, contact the customer service.
- This device must only be supplied at safety extra low voltage (SELV) corresponding to the marking on the appliance.
- To disconnect the unit from the mains, always pull the plug out of the power connection. Do not pull on the cable!

### Specific groups of people

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children must be supervised to ensure that they do not play with the device.
- Do not use this unit to supplement or replace medical treatment. Chronic complaints and symptoms may be made worse.
- Do not use the NM 990 neck massager, or only after consulting your doctor, if:
  - you are pregnant,
  - you have a cardiac pacemaker, artificial joints or electronic implants,
  - you suffer from one or more of the following illnesses and complaints: Circulation disorders, varicose veins, open wounds, bruises, broken skin, inflammation of the veins.
- Do not use the appliance near your eyes or near any other sensitive areas.
- If you experience any pain or discomfort during the massage, stop immediately and consult your doctor.
- The surface of the appliance is hot. People who are insensitive to heat must be careful when using the appliance.
- Before using the massager, consult your doctor regarding any unexplained pain, if you are undergoing medical treatment or using medical equipment.

### Before using the unit

- Check the massager carefully for damage before each use. Do not operate a defective device.
- Do not use the device if any damage is visible, if it is not working properly, if it has been dropped or if it has become wet. To avoid hazards, send the unit to the service centre for repair.

### Operating the unit

- **The device is not to be used for commercial use or medical purposes.**
- **If you have health concerns, consult your doctor before using the device.**
- **Only use the massager indoors.**
- **Do not use the massager in damp rooms (e.g. when bathing or showering).**
- The device must only be used for its intended purpose as described in the instruction manual.
- Using the device for any other purpose invalidates the warranty.
- Only use the device in the proper place, that is, around the neck and on the shoulders.
- Do not use the device while you are asleep or when you are lying in bed.
- Do not use the device before going to sleep. The neck massage has a stimulating effect
- Never use the device while driving a vehicle or operating a machine.
- The maximum operating time for one session is 15 minutes. Prolonged use reduces the service life of the device.
- Very long periods of use can lead to overheating. Always allow the device to cool down before using it again.
- Never put the unit down next to an electric stove or other sources of heat.
- Never leave the unit unattended if it is connected to the power supply.

- Do not allow the massage cushion to come into contact with pointed or sharp objects.
- Do not cover up the appliance when it is switched on. Never use it under blankets or cushions. There is a risk of fire, electric shock or injury.

#### Maintenance and cleaning

- The unit is maintenance-free. You may only carry out cleaning work on the unit yourself. In the event of malfunctions or damage, do not repair the unit yourself, as this will invalidate any warranty claims. Consult your specialist dealer and have repairs carried out only by authorised service centres in order to avoid hazards.
- Only use accessories intended for the unit and supplied by the manufacturer.
- Cleaning and user maintenance must not be carried out by children without supervision.
- Before cleaning the appliance, make sure that the appliance is switched off and the power supply unit is disconnected. Allow the unit to cool down.
- Clean the massager with a slightly damp sponge only. Never use brushes, aggressive cleaning agents, petrol, thinners or alcohol. Rub the seat dry with a clean dry cloth.
- The neck cover **1** can be removed (zipper) and washed in the washing machine with mild detergent in the gentle wash cycle at 30°C.
- May not be bleached, ironed or dry cleaned!
- Never immerse the appliance in water to clean it and ensure water cannot get into the appliance.
- Wait until the device is completely dry before using it again.
- Store the device in a safe, clean, cool and dry, well-ventilated place. Secure it against falling and shocks.
- You should ideally store the device in its original packaging and keep it in a clean, dry place.

#### Charging the battery

- Use a 5V mains adapter with an output current of 2A. Connect the supplied USB charging cable to the adapter and to the USB-C port of the device.
- The LED **1** indicates the charging process by flashing slowly. When the battery is fully charged, the LED light up constantly. The device cannot be used during the charging process. You can recharge the battery at any time.
- If the LED **1** flashes quickly the remaining charge is low – the battery should be recharged as soon as possible.

Power supply: Li-Ion battery 7.4V, 2200mAh,  
Input: USB-C, 5V= 2.0A

**Disposal:** This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Consult your municipal authority or your dealer for information about disposal.

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